Delta Instrument localizzata.....

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100 campioni/h
Manuale
Led Flow Cytometer

Analytical Science in Dairy

Mercato Combiscope Central Milk testing (CMT)

DeltaInstruments
An Advanced Instruments Company

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Patrick Cree
Hendrik de Vries

Analytical Science in Dairy

Specification Combiscope FTIR 600/300 HP/CP

DeltaInstruments
An Advanced Instruments Company

Additional test results full spectrum calibration models

Fatty acids in milk	Range	SD(y)	SECV	RPD
Saturated %mm	0.13 - 4.09%	0.65	0.031	20.9
Unsaturated %mm	0.06 - 1.98%	0.29	0.035	8.5
PolySaturated%mm	0.05 - 1.76%	0.25	0.034	7.6
MonoSaturated%mm	0.01 - 0.24%	0.04	0.017	2.5
C16:0 %mm	0.06 - 1.88%	0.28	0.042	6.6
C18:0 %mm	0.02 - 0.83%	0.14	0.043	3.3
C18:1 trans %mm	0.04 - 1.53%	0.21	0.035	6.0

SECV = standard error of cross validation, cyclic and leaving out #1 sample at the time. RPD = SD(y)/SECV

IDF and AOAC approved methods

Component	Range	Performance range	Repeatability	Accuracy bulk	Accuracy single Cow
Fat (A)	0-15%	2-15%	Cv<0.5%	Cv<1%	Cv<2.0%
Fat (B)	0.15%	2-15%	Cv<0.5%	Cv<1%	Cv<1.5%
Protein	0-10%	2-15%	Cv<0.5%	Cv<1%	Cv<1.5%
Lactose	0-10%	2-15%	Cv<0.5%	Cv<1%	Cv<1.5%
Solids	0-20%	2-20%	Cv<0.5%	Cv<1%	Cv<1.5%

Somatic Cell counting

Range 0-10000K cells/ml, performance range 0-2000K cells

Level	Repeatability	Accuracy
100K cells/ml	Cv<6%	<10%
100K-500K cells/ml	Cv<4%	<10%
500K-2000K cells/ml	Cv<3%	<10%

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
Methane (CH_4) – greenhouse gas emissions

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

Analytical Science in Dairy

Pacchetto Adulteranti

DeltaInstruments
An Advanced Instruments Company

Protein (02 01)	parm021=	Protein standard milk model
NPN/CU (05 01)	parm051=	NPN/CU in mg urea/100g,
NPN/CU (05 02)	parm052=	NPN/CU as M.U.N./100g = in mgN/100g,
NPN/CM (05 60)	parm056=	NPN/CM in mg melamine/kg,
NPN/CA(05 70)	parm057=	NPN/CA in mg Amonium NH3/kg,
Conductivity (08 02)	parm082=	C*, sample conductivity in mS, *: C = max 20mS
Citrate (11 01)	parm111=	Citric acid/Citrate in milk PLS model
FFA (12 02)	parm122=	Free Fatty Acids
pH (14 01)	parm141=	pH of raw milk valid in the pH-range: 5.4-6.6
Sucrose (61 03)	parm613=	Sucrose %m/m milk added sugar
Melamine (73 00)	parm730=	Melamine in mg/kg
Melamine pH (73 01)	parm731=	Melamine pH corrected in mg/kg derived parameter
Fat (95 07)	parm957=	Fat (Comprehensive Global Fat model)
Sat Fat (96 01)	parm961=	Saturated Fat %m/m
Unsat Fat (96 02)	parm962=	Unsaturated Fat %m/m
Cis Fatty Acids (97 09)	parm979=	Cis Fat %m/m
Trans Fatty Acids (98 00)	parm980=	Trans Fat %m/m

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

Analytical Science in Dairy

Delta Instruments
An Advanced Instruments Company

DeltaInstruments
An Advanced Instruments Company

DeltaInstruments
An Advanced Instruments Company

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- **USA Cornell University, Dave Barbano**

CombiScope FTIR 400

- **Brazil University Federal Minas Gerais,**

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
Silvia Orlandini

Somascope Smart

DeltaInstruments
An Advanced Instruments Company


Analytical Science in Dairy

Progetto IDF/ICAR Reference System on Somatic Cell

DeltaInstruments
An Advanced Instruments Company

Routine somatic cell counting in milk is largely performed with automated instruments. The accuracy of the counting level of these worldwide spread instruments is checked with various available reference materials. It is noted that most reference material providers are involved in intercomparisons with other reference material providers. These intercomparisons are very valuable as a beginning of a globally interlinked networked reference system. The implementation of such a global system for somatic cell counting in milk is the challenging target of the joint IDF/ICAR project group.*

* ICAR = International Committee for Animal Recording

Outcome of questionnaire for routine laboratories

By Silvia Orlandini*, Véronique Ninane** and Ute Braun***

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Achieving better equivalence in somatic cell counting requires a thorough understanding of the current status with applied systems for ensuring the accuracy of the analyses performed worldwide. The questionnaire dispatched to reference material providers in early 2010 brought a good overview of the available reference materials. A second questionnaire on the use of these reference materials and the calibration of routine instruments was dispatched to routine laboratories in November 2010. The main outcome of the 205 (!) received replies is summarized below.

Analytical context

As a first part, a brief inventory was made on the type of instruments and technologies, the tools for control and the calibration procedures. From the answers received it appeared that there is confusion between the terms calibration and calibration adjustment. Adjustment does not seem to be practiced extensively, only 30 % of the laboratories reported a limit of tolerance for slope and bias deviations. Many laboratories seem simply to recalibrate, even if the deviations are within tolerances.

Analytical Science in Dairy

DeltaInstruments
An Advanced Instruments Company

Delta Instruments

Analytical science in Dairy

Analytical Science in Dairy